

Outer Thigh Abductor-HM019

- | | |
|-----------------------|---|
| 1: Main Frame Gauge | "Flat Oval Tube 50*100*2.5mm; Shroud Frame Rectangular Tube 50*110*2.0mm; |
| | Motion Arms Flat Oval Tube 40*80*2.0mm" |
| 2: Wire Spec. | Dia 5.6mm |
| 3: Cushion | Material: High strength sponge, artificial PVC leather |
| 4: Guiding Bar | Round 20*2.0mm Stainless Steel Bar |
| 5: Shroud Cover | Integrated iron guard plate makes the equipment safe, beautiful and seamless |
| 6: Handle | Made of TPV material, soft, non-slip and durable |
| 7: Adjusting Position | Pulling bar adjustment structure |
| 8: Weight Stack | Steel plate thickness : 18mm; Better impact and wear resistance; Standard configuration: 80kgx1 group, 5kg/piece |
| 9: Pulleys | High Strength Nylon Pulleys, which have passed the over 100,000 times non-stop fatigue tests |
| 10: Accessories | High strength bolts and locknuts, which can meet National Standard Grade 8.8, are better in strength and rust resistance. |
| 11: Size/ N.W. | Size: 1349*1250*1601mm (L*W*H) |
| 12: Size/ N.W. | |
| 13: Manual | Chinese Manual English Manual |
| 14: Muscles Trained | Exercise the thigh abductor muscles. |